



HEMI

Modulares Armschalensystem
 Modular Arm Support System
 Modulair armschaalsysteem
 Modulært armstøttesystem
 Système modulaire d'accoudoir concave
 Sistema a coppa per braccio modulare

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01. Purpose and Indication

The Modular Arm Support System HEMI is normally used in the nursing, hospital and rehabilitation area. It is used for arm positioning (e.g. in case of hemiparesis or hemiplegia) on the wheelchair and fixed with the fastening unit U-BEF-HEMI-P. These instructions for use are intended for qualified personnel in the fields of rehabilitation and care.

02. Contraindication

The Arm Support System is unsuitable for all applications, which deviate from the intended use. The HEMI armrest system can only be used if it can be installed since part of the wheelchair. The installation may only be carried out by qualified specialist personnel.

03. Safety Instructions

This operating instructions must be read and observed before using.

- The Arm Support System is only intended for use on the side part of the wheelchair. It may only be installed and adjusted by qualified specialist personnel. Check that the armrest system is firmly seated before use.
- The operating instructions are part of the HEMI Arm Support System. It must be kept and if it is passed on, it must be handed over.
- The HEMI Arm Support System must not be used to carry, steer or push the wheelchair.
- Replace damaged parts of the armrest system immediately! If there are visible defects, the HEMI Arm Support System may no longer be used.
- Observe the requirements for cleaning and disinfection in section 5.
- Keep away from hot objects and sources of ignition. The PU foam parts could be damaged.
- If serious incidents occur in connection with this product, the incidents must be reported to Rolko Kohlgrüber GmbH and to the competent authority of the EU Member State in which the user or patient is based.

04. Assembly

The operating instructions must be read carefully before installation. Please check the items for damage. The assembly may only be carried out on wheelchairs with suitable side parts. The Arm Support System can be put together individually. Each part must be ordered separately.

1. Assembling the Arm Support System (Fig. 1): Place the forearm shell (PU-HEMI-1-S) and the desired handrest (PU-HEMI-2-S, PU-HEMI-3R-S, PU-HEMI-3L-S or PU-HEMI-5-S) on a table with the profiles facing upwards and the open side together. Insert the connecting plate halfway into the rail of the palm rest so that the two holes for the grub screws protrude. Now slide the forearm cup onto the connecting plate. Screw the 4 setscrews into the 4 holes of the connection plate. Optionally, the extension piece can be mounted between the forearm shell and the desired palm rests. When assembling, proceed exactly as just described.
2. Mounting the fastening unit (Fig. 2): The fastening unit (U-BEF-HEMI-P) is suitable for wheelchair side parts with round tube \varnothing 19 mm and flat surfaces. Hole spacings of approx. 16 mm to 140 mm are possible. Put the supplied half clamps and washers M5 / M6 together and screw both with M5 / M6 screws through the holding plate on one side part of the wheelchair. If the assembly is carried out on flat surfaces, the half clamps are not required. **Please note that the M5 / M6 screws are not included in the scope of delivery.**
3. Fixing the armrest on the fastening unit (Fig. 3a): Insert the 2 slot nuts in the rail of the forearm shell so that the holes in the slot stones and attachment unit are one above the other. Fasten the armrest to the fastening unit in the desired position using the Allen screws supplied. The rotation of the armrest is individually adjusted using locking bolts.

locking pin

05. Cleaning and Desinfecting

Clean the armrest system with lukewarm water with the addition of a neutral detergent. Do not use detergents and disinfectants that contain concentrated acids and bases, chlorinated hydrocarbons, chlorine and alcohols. Also avoid contact with hot water, water vapor and technical oils.

06. Maintenance

Check the individual parts of the armrest system carefully for optical defects and damage before use. Worn or damaged parts must be replaced immediately. Regularly check that the fastening unit is firmly seated on the wheelchair and the function of the armrest system.

07. Re-use

The product is suitable for re-use. Due to the influence of weather and UV radiation, mechanical load such as natural ageing or careless handling can cause premature abrasion on which the manufacturer has no direct influence. If the surface or the mounting unit is damaged, dispose of the product immediately. The fastening unit, PU armrests and accessories such as sliding blocks can be reordered from specialist retailers. Basic cleaning is required for re-use, see section 5: "Cleaning and Desinfecting"

08. Disposal

The fastening unit consists of aluminium, stainless steel and steel. This can be disassembled and the individual parts can be disposed of separately. The foam parts are made of polyurethane and aluminum. These can be disposed of with household waste.


09. Technical Data


Item description:

Forearm cup (incl. two sliding blocks)
 Hand rests (incl. connecting element and four grub screws)
 - flat
 - ergonomic right
 left
 universal

Item no.:	Length:	Width:	Height:
PU-HEMI-1-S	285 mm	123 mm	75 mm
PU-HEMI-2-S	180 mm	136 mm	32 mm
PU-HEMI-3R-S	192 mm	152 mm	
PU-HEMI-3L-S	192 mm	152 mm	
PU-HEMI-5-S	241 mm	166 mm	
PU-HEMI-4-S	70 mm	95 mm	32 mm
U-BEF-HEMI-P			
U-BEF-HEMI-P-2			

Extension piece (optional)
 Swivelling fixation unit (incl. two sliding blocks, two hexagon screws M5/M6 and mounting plate for wheelchair attachment)
 Alternative mounting plate with 37.5 mm outreach (optional)
 Max load capacity: 40 kg

-  Manufacturer
-  Follow the instructions
-  Production date
-  Item no.
-  Batch
-  Medical device
-  Unique device identification

As an accessory for medical products of class 1, this product has a CE mark in accordance with regulation (EU) 2017/745 (MDR) and thus complies with all German and European safety requirements, that apply to the product. The product fulfils the corresponding sections of standard specification DIN EN ISO 21856. The EU declaration of conformity is available for download at instructions.rolko.com.

Operating Instructions: Arm Support System / 2024-05-27 / Rev. 06
 Basis-UDI-DI: 426067605PUHEMI00000000YG

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Reprint, also in extracts, only with prior permission by the editor. All rights reserved. The delivered product may differ from the illustrations in this manual. All technical data are characteristic and subject to construction and manufacturing conditioned tolerances.

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© Rolko Kohlgrüber GmbH, Industriestr. 14 - 16, 33829 Borgholzhausen.
CH-REP. CH-Importeur: INNOVAhelp GmbH, Thunstr. 5a, CH-3110 Münsingen