Infusionsstäbe für Rollstühle Infusion Poles for Wheelchairs Infuusstanden voor rolstoelen Infusionsstænger til kørestole Support à perfusion pour fauteuils roulants Aste portaflebo per sedie a rotelle

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Instruction Manual

Dear Customer,

the Rolko team thanks you for the trust you have placed in us by purchasing our product. With the purchase of one of our infusion poles you have chosen a highly functional product. Our infusion poles were developed taking into account the current standards, using the latest technologies and with the thoroughness and team spirit of our employees. We are certain that we can provide the best possible support for those in need of care with our high utility value and solid quality.

01. Purpose and Indication

These instructions are aimed at specialists in the care and hospital sector. The infusion poles for wheelchairs enable gravity infusion while sitting in the wheelchair. The hooks on the bottle cross are used to hold the infusions.

02. Contraindication




The infusion poles are only suitable for common wheelchairs with a round tube diameter of 22 or 25 mm. The poles are unsuitable for all other tube diameters as well as for active wheelchairs. The poles are to be attached either to the push handle or to the frame of the wheelchair. It must be ensured that no control elements, such as brake levers, are restricted in their function. The poles must be placed in the vertical position as precisely as possible. The 2 kg load capacity of the plastic hook on the bottle cross must never be exceeded and can only be guaranteed up to a deviation of 10 ° from the vertical position of the bar. The necessity of an infusion pole and the assembly on the wheelchair are the sole responsibility of trained specialists.

03. Safety Instructions




- These instructions for use must be read and followed before using the infusion poles.
- Check the infusion wand for damage before using it for the first time.
- When assembling, make sure that the infusion poles is as vertical as possible.
- It is essential to note the maximum load capacity of 2 kg per hook.
- The patient cannot move around independently in a wheelchair while the infusion is being administered. Instructions of the medical staff must be followed!
- A damaged / defective infusion poles must never be used, but need to be replaced immediately.
- The instructions for use are part of the infusion poles and must be passed on when it is passed on.
- Serious incidents in connection with this product must be reported to Rolko Kohlgrüber GmbH, as well as to the competent authority of the EU member state in which the user is established.

04. Usage

The infusion poles for wheelchairs consist of a fastening unit for 22 or 25 mm tube diameter, depending of the model with articulated joints and a 2-armed plastic bottle holder (pic. 1). Each hook of the bottle holder can hold an infusion container weighing up to 2 kg. The model MT-IN-R22 is also height adjustable. The infusion pole must be placed in a vertical position via the fastening unit or the articulated joints. The mounting unit of the MT-IN-R22 (pic. 2) enables the mounting on a horizontal tube of the frame (pic. 3). The tube is infinitely adjustable in height.

Thanks to the linkage with the two articulated joints, the fastening unit of the model MT-IN-SCHV22 (pic. 4) can be mounted on almost any suitable tube in the frame (pic. 5 and 6). Thanks to the articulated joints, it can be adjusted individually. The linkage can be detached from the fastening unit at any time by loosening the quick release and thus dismantled, while the fastening unit remains on the wheelchair.



Instruction Manual


Thanks to the infinitely variable adjustment of the inclination of the fastening unit of the model MT-IN-V22/25 (fig. 7), this infusion pole can also be attached to almost any suitable tube (\varnothing 22/25 mm) in the frame (fig. 8 and 9). The height of the infusion pole is infinitely adjustable.

Pic 7

Pic 8

Pic 9

05. Cleaning and Disinfecting

For cleaning, the infusion poles can be wiped with a damp cloth, if necessary with the addition of a mild detergent. Use lukewarm water and, if necessary, non-abrasive, pH-value neutral, commercially available household cleaners. If necessary, the infusion poles can be cleaned with commercially available disinfectants by spray or wipe disinfection. The use of cleaning agents and disinfectants with chlorine or organic or inorganic substances with active chlorine is not permitted.

06. Maintenance

Check the function of the infusion poles before each use and make sure that it is undamaged and in one place is in perfect condition.

Before each use, check:

- the tight fit of the fastenings,
- the vertical position of the infusion poles,
- the product visually for damage.

07. Re-Use

The product is suitable for reuse. Clean the infusion pole according to the specifications under "05. Cleaning and disinfection" and carry out the maintenance according to section "06. Maintenance".

08. Disposal

All parts of the infusion poles should be safely disposed of with household waste. If possible, dismantle the infusion poles into its individual parts to enable the materials to be sorted.

09. Technical Data

Item no.: MT-IN-R22 | MT-IN-SCHV22 | MT-IN-V22/25

Delivery scope: 1 infusion pole, manual


Max. load: 2 kg per hook

Temperature: 0°C up to +50°C

Used materials: aluminium, steel, stainless steel, PA6, PP

Relative humidity: max. 80%

Label:

LOT Production, batch

REF Item no,

Manufacturer

Consider operating instructions

Production date

As an accessory for medical products of class 1, this product has a CE mark in accordance with regulation (EU) 2017/745 (MDR) and thus complies with all German and European safety requirements, that apply to the product. The product fulfils the corresponding sections of standard specification DIN EN 12182. The EU declaration of conformity is available for download at instructions.rolko.com.

Instruction manual: Infusion poles - MT-IN-R22, MT-IN-SCHV22, MT-IN-V22/25 / Stand 2022-07-12 / Rev. 04

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Reprint, also in extracts, only with prior permission by the editor. All rights reserved. Alterations for reasons of technical development reserved.
All technical data are characteristic and subject to construction and manufacturing conditioned tolerances. Errors excepted.

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© Rolko Kohlgrüber GmbH, Industriestr. 14 - 16, 33829 Borgholzhausen
CH-REP. CH-Importeur: INNOVAhelp GmbH, Thunstr. 5a, CH-3110 Münsingen