TGK

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Instruction Manual

Dear Customer,

we thank you for your confidence in our company by purchasing one of our products. We are convinced that our trapeze handle for children entirely meets the requirements of the daily nursing routines. Our trapeze handle for children was developed according to the highest medical standards, using the most advanced technologies and with the thoroughness and the team spirit of our employees. We are confident that our trapeze handle makes everyday life of nursing staff and people in need of care easier and more pleasant. Please take note of instructions listed below and consider specially the safety instructions to guarantee the proper use of our product.

01. Purpose and Intended Use

Purpose: as an accessory for sick and care beds for children and adolescents, it helps with the positioning of young patients and compensates for restrictions caused by disabilities.

Intended use: this manual is aimed at trained nursing staff, the user of the handle, and the person who put the handle into operation. Pay particular attention to the safety instructions before using the handle for the first time. All drawings and graphics are not to scale and may vary.

This trapeze handle is designed for medical beds for children and adolescents in accordance with DIN EN 50637 and is only intended for use in the hospital or rehabilitation area. The handle for children is used by means of lifting supports (bed gallows) on hospital or care beds. It hangs above the bed on the bed gallows within the patient's reach. Patients can use the handle to stand up from the lying position to the sitting position on their own, provided that residual forces remain in the patient's upper body and upper extremities. Other applications must be agreed in writing with the manufacturer beforehand. The product is to be used as a work equipment for care/nursing purposes and is subject to the regulations of the responsible employer's Liability Insurance Association (in German: Berufsgenossenschaft).

02. Contraindication

The handle is unsuitable for care and hospital beds for adults in accordance with DIN EN 60601-2-52 as well as all applications that deviate from the intended use and can lead to danger. This includes, for example, the load beyond the permissible max. load (see section *Safety Instructions and Remaining Risks* and section *Technical Data*) or operation by untrained people. The handle is not suitable as a sports device for rehabilitation exercises. Pay attention to a suitable tube diameter for the loop.

03. Safety Instructions and Remaining Risks

- These operating instructions are to be read prior to use of the handle.
- The handling of this product requires exact knowledge and observation of the instructions listed below to avoid operating mistakes and grant a failure-free operation. Keep this instruction manual for the operating personnel always accessible until the disposal of the product.
- Do not use the handle against the intended purpose and avoid improper use and operation of the functions.
- Before using the handle the user needs to familiarize himself with the functional reliability of the handle and verify the proper condition of the trapeze bar.
- Check the product for damages regularly (see section *Maintenance and Inspection*).
- Replace the handle immediately if any damage to the handle and or its components are visible.
- Observe the cleaning and disinfection specifications.
- For safety reasons, avoid the permeation of small foreign bodies and/or liquids in the handle (see section: Care, Cleaning, Disinfecting).
- Keep the handle away from strong heat sources such as cigarettes or open flames and observe the specified ambient conditions below.
- Do not exceed the safe working load of 45 kg.
- If serious incidents occur in connection with this product, the incidents must be reported to Rolko Kohlgrüber GmbH and to the competent authority of the EU Member State in which the user or patient is based.
- This handle is unsuitable for nursing and hospital beds for adults.

Instruction Manual

04. Operation

Read this instructions manual carefully and check if the product is complete or has any damages. An initial cleaning of the handle is necessary, in compliance with the permitted cleaning substances and disinfectants (see section: *Care, Cleaning and Disinfecting*), this also applies for the further uses of the product. Perform the mounting on the bed gallows as described in the following instructions.

- **Attachment to the trapeze bar**

Hang the handle with the loop on the bed gallows. The loop must be properly attached to the attachment points of the bed gallows and the anti-slip area must be between the belt strap and the bed gallows.

**The attachment is only allowed on a trapeze bar, which is intended for this handle.
Please also consider the instructions notes of the bed and bed gallows**

- **Application and use**

The handle may be used with one or both hands.

Do not overload the handle (max. 45 kg)!

- **Park position when not in use**

While the handle is not being used, you may hang it on the bed gallows. It must be ensured that the handle can not accidentally slip off the bed gallows.

Instruction Manual

05. Care, Cleaning and Disinfecting

Before cleaning and disinfecting the handle, remove it from the bed gallows / trapeze bar. Please consider the following guidelines when cleaning.

The cleaning and disinfection of the product must not be carried out in pressure chambers / autoclaves, sterilisation, cleaning or disinfection machines!

Clean and disinfect the handle before the first use and also before each reuse. Disinfect the handle on a regular basis and at least prior to each change of patient. The shape of the handle is designed to facilitate the cleaning. In addition, the handle is made of polypropylene and thus it has a high chemical resistance, which makes the hygiene process easier. There were no additional biocidal agents used during the manufacturing of the product so that allergy and resistance risks were not increased.

• Care and Cleaning

For an easy cleaning the trapeze handle may be wiped clean with a damp cloth and if necessary with the addition of a mild detergent. Use lukewarm water and if necessary, non-abrasive, pH-neutralized and commercial household cleaners. Please also note the following section *Disinfection / Sterilisation*. After each cleaning and disinfection please let it dry completely so that no residual moisture remains in the webbing or other parts.

• Disinfection / Sterilisation

We recommend the detergents and disinfectants listed below to clean the handle. The usage must be carried out only according to manufacturer specifications! Observe the corresponding product information, in particular regarding the dosage and a possible hazard in combination with other cleanings / disinfection agents. Before using other detergents or disinfectants please consult the agent manufacturers in any case. Due to the variety of detergents and disinfectants available on the market, as well as constant innovation, we can not guarantee any kind of compatibility. In case of doubt please contact the manufacturer, since a misapplication may lead to a premature wear and failure of the material presenting a danger to the user.

- From the field of formaldehydes and/or other aldehydes respective derivatives: ALDASAN 2000, B5, BUDENAT ACUTE D 441, DESIFOR PROTECT, DESOMED PERFEKT, HYGENIC VIRUZID, KOHRSOLIN EXTRA, LYSOFORM, LYSOFORMIN, MELSITT, MINUTIL, NÜSCOSEPT
- From the field of percompounds: APESIN AP100, APESIN AP 100 PLUS, DISMOZON PLUS, DISMOZON PUR, INCIDIN ACTIVE, PERFORM, TERRALIN PAA, WOFASTERIL
- From the field of carbolic acid (derivate): AMOCID

The following substances may **not** be used, since they alter the material properties, leading to premature wear, which may cause a potential material failure and thus hazards to the user:

- All phenol or phenol derivatives not mentioned in this manual, as well as other aromatic hydrocarbons.
- Chlorin organ. or anorganic substances with active Chlorin.
- Alkalies.
- All formaldehyde and / or other aldehydes or derivatives not mentioned in these instructions for use.

06. Maintenance and Inspection

The lifetime of the handle including the belt strap is influenced by many factors, such as temperature, humidity, UV radiation, exposure to chemicals, mechanical stress, etc. Also, a careless handling may cause wearout and premature aging of the plastic parts, a manufacturer has no influence in this process. Therefore, we recommend to exchange the product after every 5 years as a precaution, even if no signs of external damage exist. The recommended service life, as well as the maintenance intervals stated below are valid from the date of commissioning.

• Maintenance intervals

Simple visual and functional check by the nursing staff before each use on a new patient.
Functional test, control of the belt strap, the handle, the anti-slip part and the seams: annually.

• General test

Inspect the handle regularly regarding safe functionality and visual impairments. All visible parts must be inspected regularly for wearout or damage marks, for example, by scratches, discolorations, cracks etc. The webbing and stitching must not be cut, nor look twisted or frayed.

If the trapeze handle and/or its components are damaged or show any signs of malfunction, the product must be replaced immediately!

Instruction Manual

07. Storage and Reutilization

Before storing the handle the following steps need to be performed:

1. Conduct a function test before storing the product.
2. Cover or wrap the handle so that it does not get dirt or damaged. Protect from prolonged/permanent UV light.
3. Due to the maintenance periods attach a date of storage to the product.
4. Ensure that this operating instructions are present.
5. Store the handle only in closed rooms with a temperature between +10 and +40 °C and with a max. humidity of 80%. Before reusing it please conduct the cleaning and disinfection steps as described under section: *Care, Cleaning and Disinfecting*.

08. Disposal

The entire product must be disposed properly and professionally, in accordance with the laws and regulations of your region.

09. Technical Specifications

Dimensions of the handle (height x width): ca. 220 x 265 mm

Height adjustment range (lower edge from grip area to the upper edge of the loop): ca. 470 mm (may vary depending on the model)

Max. load capacity: 45 kg


Own weight: ca. 271 g (may vary depending on the model)

Materials used: PP, PES and TPE

Permitted ambient conditions:

Temperature between +10 °C bis +40 °C

Relative humidity up to max. 80%

LOT production batch
MD REF Medical device, item no.
manufacturer
consider operating instructions

As an accessory for medical products of class 1, this product has a CE mark in accordance with regulation (EU) 2017/745 (MDR) and thus complies with all German and European safety requirements, that apply to the product. The product fulfills the corresponding sections of standard specification DIN EN 50637. The EU declaration of conformity is available for download at instructions.rolko.com.

Operating Instructions: Medical Trapeze Handle for Children / TGK-1 / TGK-RGF33 / 2023-09-26 / Rev. 6
Basis-UDI-DI: 426067605TGK0000000000R7

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Reprint, also in extracts, only with prior permission by the editor. All rights reserved.

The delivered product may differ from the illustrations in this manual.

All technical data are characteristic and subject to construction and manufacturing conditioned tolerances.

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© Rolloko Kohlgrüber GmbH, Industriestr. 14 - 16, 33829 Borgholzhausen