

Designed Wheelchair Spoke Guards

One of the most remarkable parts of a wheelchair are the spoke wheels. Why not give your wheelchairs a colourful touch of personality with your own logo or one of our designs on a spoke guard, while keeping the fingers of the wheelchair users out of danger? Your company logo on a spoke guard also helps to identify your wheelchairs in e.g. hospitals, nursing homes or medical institutions.

- Stable and durable material, thickness: 1 mm or 1.5 mm.
- Suitable for all common spoke wheels.
- Scratch-resistant and high-quality printing on the backside of the spoke guard.
- Attachment on the wheels with cable ties or clips.
- Spoke guard colour: transparent or completely coloured, for example in white.
- You do not have an artwork? Upon request we provide you with one of our professional artworks or with your company logo. Please contact us for more information.



Enquiry form on the the backside ->

Spoke Guard Sizes



For 4/6 position push rim fixation

| Thickness (mm) | Zoll Inch | Art-Nr. Art.-No. |
|----------------|-----------|------------------|
| 1 | 20" | SP-TR-510L4 |
| 1 | 22" | SP-TR-570L4/6 |
| 1 | 24" | SP-TR-620L4/6 |
| 1,5 | 24" | SP-TR-620-1 |
| 1 | 26" | SP-TR-670L4/6 |



For 6/12 position push rim fixation

| Thickness (mm) | Zoll Inch | Art-Nr. Art.-No. |
|----------------|-----------|------------------|
| 1 | 28" | SP-TR-720L6/12 |



Without fixation notches

| Thickness (mm) | Zoll Inch | Art-Nr. Art.-No. |
|----------------|-----------|------------------|
| 1 | 20" | SP-TR-20x1.75 |
| 1 | 22" | SP-TR-570 |
| 1 | 24" | SP-TR-620 |

Enquiry

Please send your enquiry to fax number: **+49 (0) 5425 9402-21** or e-mail address: info@rolko.com
(Minimum order quantity for spoke guards with your own design: 50 pieces. There is no minimum order quantity for our designs)

| Your Design | | Quant. | Art. no.* |
|---|------------------------------|-----------------|--------------------------------|
| Transparent background | Coloured background | | _____ |
| | | | *See above for article numbers |
| Number of colours: | | Client no.: | _____ |
| Please send your logo without any commitment to: info@rolko.com | | Company: | _____ |
| | | Contact person: | _____ |
| | | Address: | _____ |
| | | | _____ |
| Rolko Design | | Phone: | _____ |
| Design name: | _____ | Fax: | _____ |
| Size: | 22" 24" | E-mail: | _____ |

Do you have questions? Then please contact us! We are happy to help.